



## iCare tonometer instructions

Written by Administrator

Wednesday, 09 October 2013 07:48 - Last Updated Wednesday, 09 October 2013 08:52

---

from the VERY CENTRE of the cornea.

4. The probe tip should be 4 to 8 mm from the cornea. This can be tricky to set up yourself, and is the reason for the review in clinic initially, to be sure that this is correct.

5. Press the measuring button gently. Remember, keep the iCare machine steady while you press the button. Practise this with your hand held up well away from your eye before trying for the first time.

6. Six measurements are taken. A longer beep is sounded once all measurements are collected, and then the average is displayed on the machine, with the IOP range given.

The tonometer turns itself off automatically after one minute of inactivity.